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



<p>実験1: 丸底フラスコ</p>	 <p>The screenshot shows the T Fab Graph app interface. On the left, there is a green sidebar with a schematic of a round-bottom flask connected to three sensors (blue, pink, yellow). Below the schematic are controls: a toggle for 'AkaDakoを切断する' (Turn off AkaDako), and checkboxes for '対応センサー' (Corresponding sensor), 'プライバシーポリシー' (Privacy Policy), and 'AkaDako.com'. The main area displays a line graph titled '14.水温(デジタルA1)[°C]' (14. Water temperature (Digital A1) [°C]). The graph shows a temperature curve that starts at approximately 20°C and rises to a plateau at 100°C. The current temperature is 99.3°C. The x-axis represents time from 16:40:20 to 16:52:36. The y-axis represents temperature from 10 to 110°C. A play button and trash icon are visible at the bottom.</p>
<p>実験2: ビーカー(蓋あり)</p>	 <p>The screenshot shows the T Fab Graph app interface for a beaker with a lid. The sidebar and controls are identical to the first experiment. The main graph shows a temperature curve rising from approximately 20°C to a plateau at 100°C. The current temperature is 100.1°C. The x-axis represents time from 17:14:04 to 17:26:40. The y-axis represents temperature from 10 to 110°C. A play button and trash icon are visible at the bottom.</p>
<p>実験3: ビーカー(蓋なし)</p>	 <p>The screenshot shows the T Fab Graph app interface for a beaker without a lid. The sidebar and controls are identical to the first experiment. The main graph shows a temperature curve rising from approximately 20°C to a plateau at 100°C. The current temperature is 99.5°C. The x-axis represents time from 17:49:17 to 18:02:49. The y-axis represents temperature from 10 to 100°C. A play button and trash icon are visible at the bottom.</p>

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